

## SL-01 Arm Sling ( Lux )



S	M	L	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm

### Features

Inner and outer surface made of sponge material which is provide extra protection supports.

### Indications

- Used for full or half dislocations in shoulder
- Used for post operative prorection and reheal
- Used to prevent frctures formation of the distall edema

## SL-01 F Arm Sling ( Mesh )



S	M	L	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm

### Features

Inner and outer surface made of perforated fabric which is air permeability and antiperspirant material.

### Indications

- Used for full or half dislocations in shoulder
- Used for post operative prorection and reheal
- Used to prevent frctures formation of the distall edema

## SL-01 B Arm Sling ( Simple )



### Features

The inner surface of the towel and upper part surface is made of spongefabric.

### Indications

- Used for full or half dislocations in shoulder

UNIVERSAL



## SL-13 Acromio-Clavicular Bandage

### Features

It is made of breathable, non-sweaty fabric. The shoulder and armpit pads made of sponge provide comfortable use.

### Indications

Acromio-clavicular dislocations and semi-dislocations

<b>S</b>	<b>M</b>	<b>L</b>
30 cm	33 cm	36 cm
33 cm	36 cm	39 cm

## SL-14 Arm Shoulder Stabilizer



UNIVERSAL

### Features

Arm sling made of sponge fabric, used to immobilize the shoulder and upper arm. With its functional design, the front velcro can be opened when necessary, allowing easy treatment and control.

### Indications

- Used for full or half dislocations in shoulder.
- Used for post operative protection and reheat.
- To provide non-surgical stabilization.
- To prevent edema formation in case of distal fracture.
- During conservative and post-operative therapy of scapula fractures.

## SL-02 Velpeau Bandage



S	M	L	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm

### Features

It's made of sponge cloth abduction for shoulder and upper arm and allows immobilization, also the bandage around the chest increase the stability

### Indications

- Used for dislocations in shoulder
- Used after surgery for protection and reheat
- Are used for Upper of the scapula and hummerus fractures rotator cuff injury

## SL-03 Arm Sling ( 30° )



UNIVERSAL

### Features

The inner surface of the towel the outer surface of the sponge cloth produced. Provides 30° abduction immobilisation shoulder. Prevent post-operative adhesions. Mobilization and rehabilitation help.

### Indications

- The impingement syndrome of post-op
- Rotator cuff tear after repair,
- Frozen shoulder under anesthesia after manipulation.
- Subocrimal bursitis and subkorakoid.
- Used after shoulder endoprosthesis.





## SL-03 D Arm Sling ( 45°- 60° )



UNIVERSAL

### Features

The inner surface of the towel the outer surface of the sponge cloth produced. Provides 45-60° abduction immobilisation shoulder. Prevent post-operative adhesions. Mobilization and rehabilitation help.

### Indications

- The impingement syndrome of post-op
- Rotator cuff tear after repair,
- Frozen shoulder under anesthesia after manipulation.
- Subacromial bursitis and subcoracoid.
- Used after shoulder endoprosthesis.



## SL-04 Clavicle Bandage



### Features

The inner part of the product is made of towel kind of material which is adjustable and provides the shoulder pulled back to the desired tension and ensures proper posture .its soft and antiperspirant the content of the material provides comfort while use.

### Indications

- Used for clavicle fracture
- Clavicle and sternoclavicular Joint Disorders

XS	S	M	L	XL	XXL
18 cm 24 cm	25 cm 30 cm	31 cm 36 cm	37 cm 42 cm	43 cm 48 cm	49 cm 54 cm



## REF-620 Shoulder Support



### Thigh Support

Made of neoprene material ,and with compression effect prevents edema.Reduces pain by reducing muscle tone .Lower-end of the strap at the velcro for easy use, there are under arm pad.

### Indications

- Semi-dislocations injuries
- Tendinitis ostararthritis symptom
- Used after prosthesis shoulder
- Used after surgery for protection

S	M	L	XL	XXL
22 cm 24 cm	24 cm 27 cm	27 cm 30 cm	30 cm 34 cm	34 cm 38 cm



## REF-621 Humerus Brace Splint



### Features

It is a shoulder orthosis made of hard sponge material and used in all situations that require immobilization of the shoulder joint in internal rotation.

### Indications

- Proximal humerus fractures
- Following proximal humeral shaft fractures
- Non-surgical, Post operative

S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm



## REF-300 Neoprene Elbow Support



### Features

Its anatomically designed from Neoprene material.

### Indications

- Protect against Trauma
- Using for compression and heat on injured or pain area.
- Protect the area during and after sport activities

S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm



## REF-301 Neoprene Elbow Support - With Epicondylitis Band



### Features

It's made of neoprene material. This anatomical fabric is designed to create compression on the muscle groups attached to the lateral epicondyle and reduce tension and pain.

### Indications

- Epicondylitis Lateral issues
- Problem with increamentation of two joint fluid.
- Using for compression on injured or pain area.

S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm

## REF-302 Neoprene Elbow Support - With Pad



### Features


Made from Neoprene material and anatomically designed. Elbow area has a special design by eva material and provides better support and protection.

### Indications

- Protect against Trauma
- Olecranon Bursitis
- Protect the area during and after sport activities

S	M	L	XL	XXL
22 cm 24 cm	24 cm 26 cm	26 cm 28 cm	28 cm 30 cm	30 cm 32 cm

## ERSA-217 Cubital Tunnel Elbow Support



### Features

It is made of neoprene material. There is a curved steel support on the inside of the elbow.

By preventing the pole joint from bending, it prevents the ulnar nerve from being stretched and relieves the symptoms.

### Indications

- Cubital tunnel syndrome
- Ulnar nerve problems

UNIVERSAL

## REF-303 Epicondylitis Band



### Thigh Support

Made of neoprene material ,and with compression effect prevents edema.Reduces pain by reducing muscle tone. Lower-end of the strap at the velcro for easy use, there are under arm pad.

### Indications

- Semi-dislocations injuries
- Tendinitis ostarartrite symptom
- Used after prosthesis shoulder
- Used after surgery for protection

UNIVERSAL



## REF-304 Epicondylitis Band



### Features

It is made of neoprene material and anatomically designed. Has a special pad with three-point silicone, which provides compression and reduce the tension in tendon. It is used for medial or lateral epicondylitis. Easy to use.

### Indications

- Semi-dislocations injuries
- Tendinitis ostarartrite symptom
- Used after prosthesis shoulder
- Used after surgery for protection

**UNIVERSAL**



## SL-906 Arm Orthosis ( Adjustable Angle Support )



### Features

This is an adjustable elbow angle orthosis manufactured anatomically in accordance with elbow shape. Its made of sponge fabric which inside is made of Terry cloths. Angle adjustment can be fulfilled easily by a pressing button. The orthosis is controlling elbow flexion and extension. Angle flexion adjustment is possible in degrees. According to the arm length, adjustment on the orthosis can be performed.

### Indications

- Helping rehabilitation in elbow contracture
- After ligament Injury
- When require to controlled the limited motion of the elbow
- Conservative treatment after elbow dislocation

S	M	L
22 cm	24 cm	26 cm
24 cm	26 cm	28 cm



# Shoulder & Arm

